

# Easy and fast Enrichment of the Aroma Components of Cheese by MMSE


## What is MMSE?

Monolithic Material Sorption Extraction (MMSE) is a novel approach for sample adsorption and extraction using monolithic hybrid adsorbent - MonoTrap. MMSE has advantages such as effective pretreatment with simple operation and high-efficiency adsorption capacity. In addition, importantly, MMSE does not require conditioning before use.


## What is MonoTrap?

MonoTrap is a state-of-the-art silica monolithic and hybrid adsorbent having a large surface area and properties based on silica, activated carbon (graphite carbon for MonoTrap TD) and Octadecyl functional group.

## Sample Preparation Procedure by MMSE-TD




Flavor components of *Blue Cheese*






- |    |                             |    |                            |
|----|-----------------------------|----|----------------------------|
| 1  | Acetaldehyde                | 16 | 2-Decanone                 |
| 2  | Butanal                     | 17 | 2-Nonanol                  |
| 3  | Ethyl Acetate               | 18 | Decanoic acid, methylester |
| 4  | Isovaleraldehyde            | 19 | 2-Undecanone               |
| 5  | 2-Pentanone                 | 20 | Butanoic acid              |
| 6  | Ethyl butyrate              | 21 | Decanoic acid, ethylester  |
| 7  | 2-Hexanone                  | 22 | Butanoic acid, 3-methyl-   |
| 8  | Isobutyl alcohol            | 23 | $\gamma$ -Caprolactone     |
| 9  | 2-Heptanone                 | 24 | 2-Undecanol                |
| 10 | Hexanoic acid, methyl ester | 25 | 2-Tridecanone              |
| 11 | Hexanoic acid, ethyl ester  | 26 | Hexanoic acid              |
| 12 | 1-Pentanol                  | 27 | 2-Pentadecanone            |
| 13 | 2-Heptanol                  | 28 | Octanoic Acid              |
| 14 | 2-Nonanone                  | 29 | <i>n</i> -Decanoic acid    |
| 15 | Octanoic acid, ethyl ester  | 30 | Dodecanoic acid            |

• Compound are identified with spectral libraries.




1	Methanethiol	13	D-Limonene	25	2-Furanmethanol
2	Ethyl Acetate	14	Acetoin	26	Acetamide
3	2-Butanone	15	Acetol	27	2-Tetradecanol
4	2-methylbutanal	16	Dimethylpyrazine	28	2-Tridecanone
5	3-methylbutanal	17	Dimethylpyrazine	29	Hexanoic acid
6	1-Propanol	18	Dimethylpyrazine	30	Dimethyl sulfone
7	Toluene	19	2-Nonanone	31	delta-Octalactone
8	Dimethyl disulfide	20	2,5-Dimethyl-3-ethylpyrazine	32	2-Pentadecanone
9	Hexanal	21	Benzaldehyde	33	Octanoic Acid
10	2-Pentenal	22	Isobutyric acid	34	delta-Decalactone
11	3-Penten-2-one	23	2-Undecanone	35	n-Decanoic acid
12	2-Heptanone	24	Butanoic acid		

• Compound are identified with spectral libraries.



Comparison of chromatograms



Area % of flavor components from each cheese

Components	Blue Cheese	Parmesan Cheese
Acids	17.7	3.73
Aldehyds	0.29	4.33
Pyrazineds	n.d.	0.24
S compounds	n.d.	0.24

**Thermal Desorption**


**System** : T-DEX 2  
**Desorb** : Temp;200 °C  
 : Time;5 min  
**Cryo** : -150deg  
**Split** : 1:0 (Desorb;1, Split;0 cc/min)  
**Injection** : 250 °C  
**Column Flow** : He, 1 mL/min

**GC/MS Conditions**

**System** : SHIMADZU GC-2010、 GCMS-QP2010Plus  
**Column** : InertCap Pure-WAX  
 0.25mmI.D. × 60m df=0.25µm  
**Column Temperature** : 40°C (5 min)→6 °C/min→250 °C  
**Detection** : MS Scan (*m/z*28.5-600)

Conditioned ready-to-use MonoTrap for TD will be delivered to you. You only have to open the ampoule bottle before sampling.

### Product Blank Data:



GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

#### **GL Sciences Inc. Japan**

22-1 Nishishinjuku 6-chome  
Shinjuku-ku, Tokyo  
163-1130, Japan

Phone: +81-3-5323-6620  
Fax: +81-3-5323-6621  
Email: [world@glsc.co.jp](mailto:world@glsc.co.jp)  
Web: [www.glsciences.com](http://www.glsciences.com)

#### **GL Sciences Inc. USA**

4733 Torrance Blvd. Suite 255  
Torrance, CA 90503  
USA

Phone: +1-310-265-4424  
Fax: +1-310-265-4425  
Email: [info@glsciencesinc.com](mailto:info@glsciencesinc.com)  
Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)

#### **GL Sciences B.V.**

Dillenburgstraat 7C  
5652AM, Eindhoven  
The Netherlands

Phone: +31-40-254-9531  
Email: [info@glsciences.eu](mailto:info@glsciences.eu)  
Web: [www.glsciences.eu](http://www.glsciences.eu)

#### **GL Sciences (Shanghai) Limited**

Tower B, Room 2003  
Far East International Plaza  
No.317 Xianxia Road, Changning District  
Shanghai, China 200051

Phone: +86-21-62782272  
Email: [contact@glsciences.com.cn](mailto:contact@glsciences.com.cn)  
Web: [www.glsciences.com.cn](http://www.glsciences.com.cn)



**International Distributors**

Visit our Website at [www.glsciences.com/distributors](http://www.glsciences.com/distributors)