


# Analysis and Retention Indices of 61 Organic Solvent Components by Nitrogen Carriers

The retention index is a relatively representative index of the retention ratio of straight-chain alkanes and is used to study constituents based on the number of carbons in the molecule. It is one of the most useful pieces of information for qualitative analysis.

The retention index can be determined because in isothermal analysis the logarithm of the retention ratio for straight-chain alkanes is linearly related to the number of carbons, and the retention ratio is also linear to the number of carbons in thermal rise analysis.

In this report, InertCap 25 was used to determine the isothermal and temperature-rise retention indices of 61 organic solvents using nitrogen as the carrier gas.

## Standard Solution Analysis



### Conditions

|             |                                                   |
|-------------|---------------------------------------------------|
| System      | : GC-FID                                          |
| Column      | : InertCap 25<br>0.25 mm I.D. × 60 m df = 0.25 μm |
| Col. Temp.  | : 40 °C - 5 °C/min - 220 °C                       |
| Carrier Gas | : N <sub>2</sub> 90 kPa                           |
| Injection   | : Split 1:50<br>240 °C                            |
| Detection   | : FID Auto Range<br>240 °C                        |
| Sample Size | : Mixed evenly<br>0.2 μL                          |

The above chromatographic conditions.

Constant-temperature analytical conditions regulate the pressure so that the linear velocity is constant.

## Retention Index in the Warming Analysis

| Peak No. | Peak name                              | Retention index | Retention time | Peak No. | Peak name                                  | Retention index | Retention time |
|----------|----------------------------------------|-----------------|----------------|----------|--------------------------------------------|-----------------|----------------|
| 1        | Methanol                               | 424             | 8.185          | 32       | 1,4-Dioxane                                | 781             | 13.158         |
| 2        | Ethanol                                | 497             | 8.498          | 33       | 4-Methyl-2-pentanone (MIBK)                | 798             | 13.598         |
| 3        | Ethyl ether                            | 522             | 8.827          |          |                                            |                 |                |
| 4        | 2-Propanol (Isopropyl alcohol)         | 529             | 8.708          | 34       | 1-Pentanol (Amyl alcohol)                  | 817             | 14.227         |
| 5        | Acetone                                | 548             | 8.930          | 35       | Isobutyl acetate                           | 821             | 14.362         |
| 6        | Tert-Butanol                           | 555             | 8.993          | 36       | Toluene                                    | 824             | 14.460         |
| 7        | Methyl acetate                         | 580             | 9.205          | 37       | 2-Hexanone (MBK)                           | 851             | 15.353         |
| 8        | Acetonitrile                           | 580             | 9.213          |          |                                            |                 |                |
| 9        | Dichloromethane                        | 595             | 9.342          | 38       | Tetrachloroethylene                        | 861             | 15.685         |
| 10       | Carbon disulfide                       | 598             | 9.368          | 39       | <i>n</i> -Butyl acetate                    | 865             | 15.847         |
| 11       | <i>n</i> -Hexane                       | 600             | 9.383          | 40       | <i>N,N</i> -Dimethylformamide              | 897             | 16.892         |
| 12       | 1-Propanol                             | 605             | 9.458          | 41       | Chlorobenzene                              | 922             | 17.802         |
| 13       | <i>Trans</i> -1,2-Dichloroethylene     | 610             | 9.552          | 42       | Ethylbenzene                               | 922             | 17.822         |
| 14       | 2-Butanol                              | 639             | 10.030         | 43       | Isopentyl acetate (Isoamyl acetate)        | 927             | 17.977         |
| 15       | Methyl ethyl ketone                    | 654             | 10.285         |          |                                            |                 |                |
| 16       | Ethyl acetate                          | 661             | 10.390         | 44       | <i>p</i> -Xylene                           | 928             | 18.032         |
| 17       | <i>Cis</i> -1,2-Dichloroethylene       | 664             | 10.437         | 45       | <i>m</i> -Xylene                           | 929             | 18.060         |
| 18       | 2-Methyl-1-propanol (Isobutyl alcohol) | 667             | 10.497         | 46       | <i>o</i> -Xylene                           | 960             | 19.188         |
|          |                                        |                 |                | 47       | Cyclohexanol                               | 961             | 19.255         |
| 19       | Chloroform                             | 675             | 10.628         | 48       | 1-Methylcyclohexanol                       | 964             | 19.338         |
| 20       | Tetrahydrofuran                        | 686             | 10.800         | 49       | Styrene                                    | 965             | 19.382         |
| 21       | 1,1,1-Trichloroethane                  | 694             | 10.937         | 50       | <i>n</i> -Pentyl acetate                   | 967             | 19.442         |
| 22       | 2-Methoxyethanol (Methyl cellosolve)   | 698             | 11.010         | 51       | 2-Butoxyethanol (Butyl cellosolve)         | 967             | 19.473         |
| 23       | Carbon tetrachloride                   | 706             | 11.202         | 52       | 2-Ethoxyethyl acetate (Cellosolve acetate) | 982             | 19.990         |
| 24       | Isopropyl acetate                      | 707             | 11.225         |          |                                            |                 |                |
| 25       | 1-Butanol                              | 711             | 11.330         | 53       | <i>N,N</i> -Dimethylacetamide              | 982             | 20.015         |
| 26       | Benzene                                | 717             | 11.488         | 54       | Cyclohexanone                              | 1000            | 20.663         |
| 27       | 1,2-Dichloroethane                     | 723             | 11.643         | 55       | 1,1,2,2-Tetrachloroethane                  | 1011            | 21.077         |
| 28       | Trichloroethylene                      | 749             | 12.318         | 56       | 4-Methylcyclohexanone                      | 1062            | 22.915         |
| 29       | <i>n</i> -Propyl acetate               | 763             | 12.683         | 57       | Phenol                                     | 1069            | 23.173         |
| 30       | 2-Ethoxyethanol (Cellosolve)           | 770             | 12.872         | 58       | 1,2-Dichlorobenzene                        | 1119            | 25.637         |
|          |                                        |                 |                | 59       | <i>o</i> -Cresol                           | 1126            | 26.082         |
| 31       | 3-Methyl-1-butanol (Isoamyl alcohol)   | 781             | 13.148         | 60       | <i>p</i> -Cresol                           | 1136            | 26.818         |
|          |                                        |                 |                | 61       | <i>m</i> -Cresol                           | 1137            | 26.868         |

\* Retention time in minutes

\* Components with a symmetry factor of 1.5 or higher are in red.


In the case of warming analysis...

Since the retention ratio of linear alkanes and the carbon number of linear alkanes are linearly related, the retention index is given by the following equation:

$$\text{Retention index } I = 100 \times \frac{t_R - t_{R(Z)}}{t_{R(Z+1)} - t_{R(Z)}} + 100 \times Z$$

 $t_R$  = retention-time of the components of interest $t_{R(Z)}$  = retention time of straight-chain alkanes emerging in front of the component of interest $t_{R(Z+1)}$  = retention time of linear alkanes emerging after the components of interestThe numbers of carbons in linear alkanes with  $Z$  = retention-time  $t_{R(Z)}$ 

Linear alkane (Z)



## Retention Index-1 in Homeothermal Analysis

| During warming<br>Peak order | Peak Name                                 | 40 °C           |                | 80 °C           |                | 120 °C          |                | 160 °C          |                |
|------------------------------|-------------------------------------------|-----------------|----------------|-----------------|----------------|-----------------|----------------|-----------------|----------------|
|                              |                                           | Retention index | Retention time | Retention index | Retention time | Retention index | Retention time | Retention index | Retention time |
| 1                            | Methanol                                  | 434             | 7.927          | 437             | 7.265          | 441             | 6.867          | 488             | 6.595          |
| 2                            | Ethanol                                   | 501             | 8.458          | 496             | 7.425          | 497             | 6.933          | 497             | 6.602          |
| 3                            | Ethyl ether                               | 532             | 8.827          | 529             | 7.533          | 526             | 6.972          | 522             | 6.617          |
| 4                            | 2-Propanol<br>(Isopropyl alcohol)         | 541             | 8.953          | 536             | 7.562          | 531             | 6.978          | 536             | 6.625          |
| 5                            | Acetone                                   | 558             | 9.235          | 557             | 7.652          | 557             | 7.020          | 559             | 6.640          |
| 6                            | <i>Tert-Butanol</i>                       | 566             | 9.380          | 562             | 7.673          | 562             | 7.027          | 559             | 6.640          |
| 7                            | Methyl acetate                            | 587             | 9.805          | 583             | 7.783          | 580             | 7.060          | 583             | 6.657          |
| 8                            | Acetonitrile                              | 586             | 9.785          | 585             | 7.793          | 589             | 7.078          | 591             | 6.663          |
| 9                            | Dichloromethane                           | 595             | 10.002         | 599             | 7.875          | 605             | 7.113          | 614             | 6.683          |
| 10                           | Carbon disulfide                          | 590             | 9.882          | 607             | 7.930          | 624             | 7.157          | 650             | 6.720          |
| 11                           | <i>n-Hexane</i>                           | 600             | 10.128         | 600             | 7.883          | 600             | 7.102          | 600             | 6.670          |
| 12                           | 1-Propanol                                | 608             | 10.352         | 604             | 7.910          | 603             | 7.108          | 608             | 6.677          |
| 13                           | <i>Trans-1,2-Dichloroethylene</i>         | 611             | 10.418         | 615             | 7.982          | 620             | 7.147          | 626             | 6.695          |
| 14                           | 2-Butanol                                 | 648             | 11.683         | 644             | 8.210          | 642             | 7.202          | 647             | 6.717          |
| 15                           | Methyl ethyl ketone                       | 662             | 12.303         | 660             | 8.358          | 660             | 7.252          | 662             | 6.735          |
| 16                           | Ethyl acetate                             | 670             | 12.690         | 664             | 8.397          | 660             | 7.252          | 656             | 6.728          |
| 17                           | <i>Cis-1,2-Dichloroethylene</i>           | 666             | 12.498         | 672             | 8.480          | 680             | 7.317          | 687             | 6.767          |
| 18                           | 2-Methyl-1-propanol<br>(Isobutyl alcohol) | 675             | 12.913         | 671             | 8.468          | 668             | 7.278          | 672             | 6.747          |
| 19                           | Chloroform                                | 677             | 13.032         | 683             | 8.597          | 689             | 7.348          | 699             | 6.783          |
| 20                           | Tetrahydrofuran                           | 685             | 13.477         | 691             | 8.695          | 699             | 7.385          | 707             | 6.795          |
| 21                           | 1,1,1-Trichloroethane                     | 691             | 13.800         | 698             | 8.792          | 709             | 7.423          | 718             | 6.812          |
| 22                           | 2-Methoxyethanol<br>(Methyl cellosolve)   | 696             | 14.158         | 700             | 8.807          | 708             | 7.418          | 709             | 6.798          |
| 23                           | Carbon tetrachloride                      | 701             | 14.450         | 712             | 8.980          | 725             | 7.490          | 737             | 6.842          |
| 24                           | Isopropyl acetate                         | 713             | 15.290         | 706             | 8.892          | 700             | 7.388          | 693             | 6.775          |
| 25                           | 1-Butanol                                 | 716             | 15.492         | 713             | 8.993          | 712             | 7.432          | 710             | 6.800          |
| 26                           | Benzene                                   | 714             | 15.343         | 724             | 9.158          | 736             | 7.540          | 749             | 6.862          |
| 27                           | 1,2-Dichloroethane                        | 723             | 16.035         | 729             | 9.237          | 737             | 7.545          | 746             | 6.858          |
| 28                           | Trichloroethylene                         | 747             | 18.265         | 756             | 9.730          | 767             | 7.697          | 778             | 6.917          |
| 29                           | <i>n-Propyl acetate</i>                   | 771             | 20.922         | 765             | 9.912          | 762             | 7.667          | 757             | 6.877          |
| 30                           | 2-Ethoxyethanol<br>(Cellosolve)           | 772             | 21.102         | 774             | 10.107         | 782             | 7.785          | 782             | 6.927          |
| 31                           | 3-Methyl-1-butanol<br>(Isoamyl alcohol)   | 785             | 22.872         | 783             | 10.320         | 789             | 7.828          | 783             | 6.928          |

\* Retention time in minutes

\* Components with a symmetry factor of 1.5 or higher are in red.

In the case of homeothermal analysis....

Since the logarithm of the retention ratio of the linear alkanes and the number of carbons of the linear alkanes are in a linear relationship, the retention index is given by the following equation:

$$\text{Retention index } I = 100 \times \frac{\log t'_R - \log t'_{R(Z)}}{\log t'_{R(Z+1)} - \log t'_{R(Z)}} + 100 \times Z$$

tR = retention time of the components of interest.

tR(Z) = retention time of straight-chain alkanes emerging in front of the component of interest

tR(Z+1) = retention time of linear alkanes emerging after the components of interest

The numbers of carbons in linear alkanes with Z = retention-time tR(Z)

t'R = corrected retention-time t'R = tR - t0

t0 = hold-up time (elution time of non-retained components)

## Retention Index-2 in Homeothermal Analysis

| During warming<br>Peak order | Peak Name                                  | 40 °C           |                | 80 °C           |                | 120 °C          |                | 160 °C          |                |
|------------------------------|--------------------------------------------|-----------------|----------------|-----------------|----------------|-----------------|----------------|-----------------|----------------|
|                              |                                            | Retention index | Retention time | Retention index | Retention time | Retention index | Retention time | Retention index | Retention time |
| 32                           | 1,4-Dioxane                                | 780             | 22.113         | 785             | 10.362         | 795             | 7.867          | 802             | 6.970          |
| 33                           | 4-Methyl-2-pentanone (MIBK)                | 798             | 24.905         | 798             | 10.708         | 802             | 7.908          | 801             | 6.967          |
| 34                           | 1-Pentanol (Amyl alcohol)                  | 821             | 28.912         | 818             | 11.275         | 837             | 8.057          | 818             | 7.008          |
| 35                           | Isobutyl acetate                           | 872             | 41.963         | 822             | 11.402         | 837             | 8.057          | 814             | 6.997          |
| 36                           | Toluene                                    | 816             | 28.002         | 828             | 11.590         | 869             | 8.215          | 855             | 7.107          |
| 37                           | 2-Hexanone(MBK)                            | 853             | 36.477         | 853             | 12.515         | 896             | 8.357          | 855             | 7.107          |
| 38                           | Tetrachloroethylene                        | 848             | 35.273         | 864             | 12.970         | 915             | 8.603          | 895             | 7.240          |
| 39                           | <i>n</i> -Butyl acetate                    | 812             | 27.313         | 867             | 13.108         | 905             | 8.445          | 858             | 7.118          |
| 40                           | <i>N,N</i> -Dimethylformamide              | 893             | 49.567         | 897             | 14.610         | 931             | 8.865          | 910             | 7.293          |
| 41                           | Chlorobenzene                              | 904             | 54.200         | 921             | 16.053         | 957             | 9.377          | 958             | 7.498          |
| 42                           | Ethylbenzene                               | 909             | 56.165         | 922             | 16.113         | 955             | 9.322          | 950             | 7.460          |
| 43                           | Isopentyl acetate (Isoamyl acetate)        | 931             | 67.440         | 927             | 16.508         | 948             | 9.177          | 920             | 7.335          |
| 44                           | <i>p</i> -Xylene                           | 914             | 58.833         | 927             | 16.492         | 959             | 9.408          | 954             | 7.482          |
| 45                           | <i>m</i> -Xylene                           | 915             | 59.048         | 928             | 16.540         | 959             | 9.422          | 955             | 7.483          |
| 46                           | <i>o</i> -Xylene                           | 942             | 73.517         | 957             | 18.787         | 982             | 9.962          | 990             | 7.663          |
| 47                           | Cyclohexanol                               | 950             | 78.497         | 959             | 18.977         | 981             | 9.940          | 987             | 7.645          |
| 48                           | 1-Methylcyclohexanol                       | 950             | 79.087         | 961             | 19.175         | 983             | 10.003         | 990             | 7.663          |
| 49                           | Styrene                                    | 948             | 77.608         | 962             | 19.253         | 984             | 10.033         | 993             | 7.680          |
| 50                           | <i>n</i> -Pentyl acetate                   | 973             | 95.490         | 968             | 19.788         | 976             | 9.807          | 960             | 7.507          |
| 51                           | 2-Butoxyethanol (Butyl cellosolve)         | 963             | 87.628         | 966             | 19.647         | 981             | 9.948          | 977             | 7.593          |
| 52                           | 2-Ethoxyethyl acetate (Cellosolve acetate) | 995             | 115.833        | 984             | 21.383         | 983             | 10.003         | 974             | 7.578          |
| 53                           | <i>N,N</i> -Dimethylacetamide              | 977             | 99.103         | 981             | 21.067         | 992             | 10.245         | 1002            | 7.730          |
| 54                           | Cyclohexanone                              | 980             | 101.488        | 994             | 22.508         | 1012            | 10.770         | 1036            | 7.945          |
| 55                           | 1,1,2,2-Tetrachloroethane                  | 999             | 119.040        | 1007            | 23.987         | 1020            | 10.947         | 1038            | 7.963          |
| 56                           | 4-Methylcyclohexanone                      | 1037            | 166.425        | 1052            | 30.552         | 1070            | 12.380         | 1092            | 8.392          |
| 57                           | Phenol                                     | 1061            | 206.025        | 1063            | 32.628         | 1072            | 12.438         | 1087            | 8.345          |
| 58                           | 1,2-Dichlorobenzene                        | 1092            | 271.508        | 1108            | 44.442         | 1121            | 15.310         | 1136            | 9.248          |
| 59                           | <i>o</i> -Cresol                           | -               | -              | 1120            | 51.382         | 1125            | 15.738         | 1134            | 9.195          |
| 60                           | <i>p</i> -Cresol                           | -               | -              | 1130            | 58.872         | 1135            | 16.843         | 1143            | 9.440          |
| 61                           | <i>m</i> -Cresol                           | -               | -              | 1131            | 59.400         | 1136            | 16.923         | 1144            | 9.460          |

\* Retention time in minutes

\* Components with a symmetry factor of 1.5 or higher are in red.

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

**GL Sciences Inc. Japan**

22-1 Nishishinjuku 6-chome  
Shinjuku-ku, Tokyo  
163-1130, Japan

Phone: +81-3-5323-6620  
Fax: +81-3-5323-6621  
Email: [world@glsc.co.jp](mailto:world@glsc.co.jp)  
Web: [www.glsciences.com](http://www.glsciences.com)

**GL Sciences Inc. USA**

4733 Torrance Blvd. Suite 255  
Torrance, CA 90503  
USA

Phone: +1-310-265-4424  
Fax: +1-310-265-4425  
Email: [info@glsciencesinc.com](mailto:info@glsciencesinc.com)  
Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)

**GL Sciences B.V.**

Dillenburgstraat 7C  
5652AM, Eindhoven  
The Netherlands

Phone: +31-40-254-9531  
Email: [info@glsciences.eu](mailto:info@glsciences.eu)  
Web: [www.glsciences.eu](http://www.glsciences.eu)

**GL Sciences (Shanghai) Limited**

Tower B, Room 2003  
Far East International Plaza  
No.317 Xianxia Road, Changning District  
Shanghai, China 200051

Phone: +86-21-62782272  
Email: [contact@glsciences.com.cn](mailto:contact@glsciences.com.cn)  
Web: [www.glsciences.com.cn](http://www.glsciences.com.cn)

**International Distributors**

Visit our Website at [www.glsciences.com/distributors](http://www.glsciences.com/distributors)