


Purity Test of Calcium Polystyrene Sulfonate, with reference to the Japanese Pharmacopoeia, 16th Edition - Comparisons of Chromosorb Series and Uniport Series of Siliceous Soil Support


Calcium polystyrene sulfonate is used as a cardiovascular drug. It reduces potassium in the blood and is used to treat hyperkalemia caused by kidney failure. In this study, a packed column was used for analysis with reference to purity test (4) described in the Japanese Pharmacopoeia, 16th edition.

Chromosorb series and Uniport series siliceous earth carriers were evaluated and compared. This report demonstrates that both supports satisfy the requirements of the test method.


Example: Measurement of standard



Example: Measurement of performance test standard of the system



Blank



Measurement conditions

Conditions

System	: GC - FID	Injection	: Direct
Column	: PEG 20M 15 % Chromosorb WAW DMCS 80/100 Uniport HP 80/100 SUS 2 m x 3 mm I.D.	Injection Vol.	: 5 µL
Col. Temp.	: 90 °C	Detection	: FID Auto Range 250 °C
Carrier Gas	: N ₂ 180 kPa	Sample	: Standard Analyte in Acetone

Measurement

Requirements for System Suitability

1. When the test is made with the system performance test standard, the number of theoretical plates and symmetry factor of the styrene peak are not less than 800 and 0.8 to 1.2, respectively.
2. When the test is repeated 6 times with the test standard, the relative standard deviation of the peak height of styrene is not more than 5%.

**Results**

Table 1. System performance test results

	Chromosorb WAW DMCS	Uniport HP
Number of theoretical plates	3502	3715
Symmetry factor	1.05	1.05

Table 2. System repeatability test results

	Chromosorb WAW DMCS	Uniport HP
1	1455	1419
2	1409	1400
3	1439	1473
4	1360	1441
5	1398	1361
6	1450	1397
Mean	1419	1415
Standard deviation	36.5	38.8
Relative standard deviation (%)	2.58	2.74

Introduction of diatomaceous earth carriers

The Japanese Pharmacopoeia contains several tests: “Siliceous earth packings for gas chromatography with XXX% liquid phase.” The two most commonly used diatomaceous earth carriers that were used in these tests are described here.

Uniport series

Uniport B (sintered after washing + acid-treated product)

Special treatment removes metals from the surface of the support, making it easy to use without catalytic activity.

It can be used with a wide range of liquid phases as well as polyglycols and polyesters.

Uniport HP (sintering, acid treatment, and silane treatment after washing)

Uniport B is specifically silanized. It is the most inert carrier for silanization.

It has good compatibility with silicon systems and is suitable to the sample that is easy to adsorb and disassemble.

Chromosorb series

Chromosorb WAW (sintered after washing + acid-treated product)

Chromosorb WAW DMCS (sintered after washing + acid treatment + silane treatment)

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
Shinjuku-ku, Tokyo
163-1130, Japan

Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
USA

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
5652AM, Eindhoven
The Netherlands

Phone: +31-40-254-9531
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
Far East International Plaza
No.317 Xianxia Road, Changning District
Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors