


Fermented foods like Miso, Sake, Syoyu made from cereals contain amino acids, organic acids, sugar and so on. In this study, they were extracted, derivatized, and analysed with gas chromatography/mass spectrometry. The silylated compounds in Miso, Sake and Syoyu extracts were identified using InertCap 5MS/NP metabolites Library (No.81).

## Procedure of sample preparation





Analysis of Miso

**GC/MS Conditions**

System : GC-MS  
 Column : InertCap 5MS/NP 0.25 mm I.D. × 30 m df = 0.25 μm  
 (Cat.No. 1010-18642)  
 Injection : Split Constant Pressure 75 kPa  
 Injection Temp. : 230 °C  
 Oven Temp. : 80 °C (2 min) - 15 °C /min - 330 °C (13 min)  
 Carrier gas : Helium  
 Septum purge flow : 5 mL/min  
 Split ratio : 25:1  
 Interface Temp. : 250 °C  
 Ion source Temp. : 200 °C  
 Ionization Voltage : 70 eV  
 Scan range : *m/z*85 - 500  
 Inj. Vol. : 1 μL

The silylated compounds in Miso extract were identified based on retention index and mass spectrum using InertCap metabolites library (Technical note No.81) and NIST.





Analysis of Sake

**GC/MS Conditions**

System : GC-MS  
 Column : InertCap 5MS/NP 0.25 mm I.D. × 30 m df = 0.25 μm (Cat.No. 1010-18642)  
 Injection : Split Constant Pressure 75 kPa  
 Inj. Temp. : 230 °C  
 Oven Temp. : 80 °C (2 min) - 15 °C /min - 330 °C (13 min)  
 Carrier gas : Helium  
 Septum purge flow : 5 mL/min  
 Split ratio : 25:1  
 Interface Temp. : 250 °C  
 Ion source Temp. : 200 °C  
 Ionization Voltage : 70 eV  
 Scan range : m/z 85 - 500  
 Inj. Vol. : 1 μL

The silylated compounds in Sake extract were identified based on retention index and mass spectrum using InertCap metabolites library (Technical note No.81) and NIST.





Analysis of Syoyu

**GC/MS Conditions**

System : GC-MS  
 Column : InertCap 5MS/NP 0.25 mm I.D. × 30 m df = 0.25 μm (Cat.No. 1010-18642)  
 Injection : Split Constant Pressure 75 kPa  
 Inj. Temp. : 230 °C  
 Oven Temp. : 80 °C (2 min) - 15 °C /min - 330 °C (13 min)  
 Carrier gas : Helium  
 Septum purge flow : 5 mL/min  
 Split ratio : 25:1  
 Interface Temp. : 250°C  
 Ion source Temp. : 200°C  
 Ionization Voltage : 70 eV  
 Scan range : m/z 85 - 500  
 Inj. Vol. : 1 μL

The silylated compounds in Syoyu extract were identified based on retention index and mass spectrum using InertCap metabolites library (Technical note No.81) and NIST.



Aloutput is available for peak identification. In that case, it is convenient to use the *n*-alkane standard solution C9-C40 in *n*-Hexane.

### What is Aloutput?

- Software for metabolome analyses based on Excel (upto 2007)
- Peak identification based on retention indices and mass spectrum
- Auto peak identification, PCA analysis, PLS analysis, t-test
- Free software!
- Everybody can download, now!

[http://prime.psc.riken.jp/Metabolomics\\_Software/Aloutput/index.html](http://prime.psc.riken.jp/Metabolomics_Software/Aloutput/index.html)

### *n*-Alkane standard solution C9-C40 in *n*-Hexane

Concentration: 50 µg/mL each component of *n*-Hexane


※C10,C20,C30,C40 are 100 µg/mL

Volume: 1 mL

Packaging: ampoule of 1 mL

#### Conditions

System : GC-FID  
 Column : InertCap 5MS/NP  
           0.25 mm I.D. × 30 m, df = 0.25 µm  
 Injection : Split 120 kPa  
 Inj. Temp. : 230 °C  
 Oven Temp. : 100 °C (2 min) - 20 °C/min - 330 °C (17 min hold)  
 Carrier gas : Helium  
 Split flow : 25 mL/min  
 Det. Temp. : 330 °C  
 Inj. Vol. : 1 µL



Description	Volume	Cat.No.
<i>n</i> -Alkane standard solution C9-C40 in C6 (cold strage)	1 mL	1021-58321
	1 mL × 5	1021-58325

This product contains hazardous material which requires special freight handling, and then additional charges apply. Please contact your local distributor.

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

#### GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome  
 Shinjuku-ku, Tokyo  
 163-1130, Japan

Phone: +81-3-5323-6620  
 Fax: +81-3-5323-6621  
 Email: [world@glsc.co.jp](mailto:world@glsc.co.jp)  
 Web: [www.glsciences.com](http://www.glsciences.com)

#### GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255  
 Torrance, CA 90503  
 USA

Phone: +1-310-265-4424  
 Fax: +1-310-265-4425  
 Email: [info@glsciencesinc.com](mailto:info@glsciencesinc.com)  
 Web: [www.glsciencesinc.com](http://www.glsciencesinc.com)

#### GL Sciences B.V.

Dillenburgstraat 7C  
 5652AM, Eindhoven  
 The Netherlands

Phone: +31-40-254-9531  
 Email: [info@glsciences.eu](mailto:info@glsciences.eu)  
 Web: [www.glsciences.eu](http://www.glsciences.eu)

#### GL Sciences (Shanghai) Limited

Tower B, Room 2003  
 Far East International Plaza  
 No.317 Xianxia Road, Changning District  
 Shanghai, China 200051

Phone: +86-21-62782272  
 Email: [contact@glsciences.com.cn](mailto:contact@glsciences.com.cn)  
 Web: [www.glsciences.com.cn](http://www.glsciences.com.cn)



International Distributors

Visit our Website at [www.glsciences.com/distributors](http://www.glsciences.com/distributors)