


InertSearch for GC

InertCap Applications

Residual solvents in pharmaceuticals (JP18 Standard stock solution Class 2 A)

Data No. GA ***-0155



Conditions

System : GC4000 Plus-FID
Column : InertCap WAX
0.32 mm I.D. × 30 m df = 0.25 µm
Col. Cat. No. : 1010-67242
Col. Temp. : 50 °C(20 min) – 6 °C/min
- 165 °C(20 min)
Carrier Gas : He 2.20 mL/min (35cm/sec)
Injection : Split 20:1
250 °C
Detection : FID Auto Range
250 °C
Sample Size : 1.0 µL
Analyte in Dimethyl sulfoxide

Analyte (R.T. min)

1. Cyclohexane	(1.619)
2. Methylcyclohexane	(1.689)
3. Tetrahydrofuran	(1.940)
4. <i>trans</i> -1,2-Dichloroethylene	(1.940)
5. Methanol	(2.081)
6. Dichloromethane	(2.238)
7. <i>cis</i> -1,2-Dichloroethylene	(2.694)
8. Acetonitrile	(2.886)
9. Toluene	(3.242)
10. 1,4-Dioxane	(3.601)
11. Ethylbenzene	(4.773)
12. <i>p</i> -Xylene	(4.982)
13. <i>m</i> -Xylene	(5.140)
14. <i>o</i> -Xylene	(6.467)
15. Chlorobenzene	(7.654)