

Simple Concentration Analysis of Lavender Flavor Components -
Using MonoTrap RCC18 (solvent-extraction)

MonoTrap RCC18 was used to simply screen and analyze volatile compounds in lavender petals. MonoTrap RCC18 is used as a collector of volatile compounds. By using MonoTrap to collect volatile compounds followed by solvent extraction, we were able to simplify the process from sampling to analysis. Using MonoTrap it was possible to detect a number of flavor components, including Linalool and Linalyl acetate in lavender.

Preliminary processing procedure

Lavender

Place 2 g of lavender petals
in a 40 mL vialCollection (HS)
One MonoTrap RCC18

60 °C for 1 hour

Solvent extraction


200 µL of
dichloromethane
5 minutes of
ultrasonication

GC/MS



GC/MS Conditions

System : GC/MS
 Column : InertCap Pure-WAX
 0.25 mm I.D. x 60 m, df = 0.5 µm
 Col. Cat. No. : 1010 - 68164
 Col.Temp. : 40 °C (5 min) - 10 °C/min - 250 °C
 Carrier Gas : He, 1 mL/min (constant flow)
 GC Inlet : 250 °C Splitless (1min)
 Detection : MS Scan (m/z 30-350)
 Sample Size : 1µL



※Standard samples are not used for qualitative analysis.
 Results from a library search.

- | | | | |
|-------------------|-----------------------|---------------------|--------------------|
| 1. Pinene | 10. Cineole | 19. Octenyl acetate | 28. Neryl acetate |
| 2. Camphene | 11. Ocimene | 20. Allo-Ocimene | 29. Bornyl acetate |
| 3. Pinene | 12. Ocimene | 21. Undecatriene | 30. Terpineneol |
| 4. Sabinene | 13. Octanone | 22. Hexyl butyrate | 31. Caryophyllene |
| 5. Carene | 14. Hexyl acetate | 23. 1-Octen-3-ol | 32. Farnesene |
| 6. Myrcene | 15. Cymene | 24. Linalool Oxide | 33. Lavandulol |
| 7. Limonene | 16. Cymene | 25. Linalool | 34. Borneol |
| 8. Butyl butyrate | 17. Carene | 26. Linalyl acetate | |
| 9. Phellandrene | 18. Hexyl isobutyrate | 27. Camphor | |

Product used

MonoTrap RCC18



Cat.No. :1050-72201

MonoTrap and Trial Kit for SE (for extracting solvents)



200 μ L glass insert (Flat Bottom)



Cat.No. :1030-17211

Image of use

InertCap Pure-WAX



Size :0.25 mm I.D. $\times$ 60 m, df = 0.5 μ m

Cat.No. : 1010-67164

* Used to elute from a rod-type MonoTrap.

- ① Set a glass insert in an autosampler vial.
- ② Place the extraction solvent in the glass insert.
- ③ Place the collected MonoTrap in the glass insert.
- ④ Add ultrapure water to a level higher than the MonoTrap in the autosampler vial
Enter the degree.
- ⑤ Irradiate the ultrasound (5 minutes).

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly.
We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
Shinjuku-ku, Tokyo
163-1130, Japan

Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
USA

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
5652AM, Eindhoven
The Netherlands

Phone: +31-40-254-9531
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
Far East International Plaza
No.317 Xianxia Road, Changning District
Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors