


## Analysis of Food synthetic sweetener

Data No. LA717-0848



### Conditions

**System** : GL-7400 HPLC system  
**Column** : Inertsil ODS-4 (5μm, 150 x 4.6 mm I.D.)  
**Column Cat. No.** : 5020-03945  
**Eluent** : 20.3 g of 10% *tetra-n*-propylammonium hydroxyde aqueous solution was dissolved in methanol:water = 20:80 (approx. 900 mL), and H<sub>3</sub>PO<sub>4</sub> was added to the solution to adjust the pH value to 4.0. Methanol: water = 20:80 was added again to make up the solution to 1000 mL.  
**Flow Rate** : 1.0 mL/min  
**Col. Temp.** : 40 °C  
**Detection** : UV 210 nm (GL-7452A PDA Detector)  
**Injection Vol.** : 20 μL

### Analyte

1. Acesulfame potassium (AK) (50 mg/L)  
 2. Sodium Saccharin (SA) (50 mg/L)  
 3. Aspartame (APM) (50 mg/L)