


Concentration Analysis of VOCs Derived from Expanded Foam Plastic Cushioning Material - Using MonoTrap a Simple Enrichment Tool

MonoTrap DCC18, a simple concentration tool was used to screen and analyze the components that volatilized from expanded foam plastic cushioning material used in a carrier case. The atmosphere in the case was collected using MonoTrap, then extracted with dichloromethane and analyzed by GC/MS. A number of peaks, including alkylbenzene compounds were detected.

Measurement procedure

Carrier case

Collection (HS)
5 MonoTrap DCC18At room
temperature
for 15 hours

Solvent extraction


500 μ L of Dichloromethane
5 minutes of ultrasonication

GC/MS

GC/MS Conditions


System	: GC - MS
Column	: InertCap Pure-WAX 0.25 mm I.D. X 30 m, df = 0.25 μ m
Col.Cat. No.	: 1010-68162
Col.Temp.	: 40 °C (5 min) - 4 °C/min - 250 °C
Carrier Gas	: He, 50 kPa
GC Inlet	: 250 °C, Splitless
Sample Size	: 1 μ L
Detection	: MS Scan (m/z 45-450)

GC/MS chromatogram



See the following page for the qualitative results of the components obtained by library searching...

Results of library search



- | | |
|--|--|
| 1. Tridecane, 4-cyclohexyl- | 26. Benzene, 1-methyl-2-(1-methylethyl)- |
| 2. 1-Propanol, 2-methyl- | 27. Benzene, 1-methyl-2-propyl- |
| 3. Propane, 1-(1,1-dimethylethoxy)-2-methyl- | 28. 2-Ethoxyethyl acetate |
| 4. Undecane | 29. Benzene, 1-ethyl-3,5-dimethyl- |
| 5. Ethylbenzene | 30. Benzene, 1,2,3-trimethyl- |
| 6. Xylene | 31. Benzene, 1-methyl-2-(1-methylethyl)- |
| 7. Xylene | 32. Benzene, 1-methyl-3-(1-methylethyl)- |
| 8. 1-Butanol | 33. Benzene, 2-ethyl-1,4-dimethyl- |
| 9. 1-Undecene | 34. Acetic acid, 2-ethylhexyl ester |
| 10. 2,3-Dimethyldecane | 35. Ethanol, 2-butoxy- |
| 11. Xylene | 36. Benzene, 1,2,4,5-tetramethyl- |
| 12. 4-Heptanone, 2,6-dimethyl- | 37. Benzene, (1-methoxy-1-methylethyl)- |
| 13. Undecane, 4-cyclohexyl- | 38. Benzene, 1,2,3,5-tetramethyl- |
| 14. Trans-Decalin, 2-methyl- | 39. 2-Ethylhexyl chloroformate |
| 15. Benzene, propyl- | 40. 2-Propenoic acid, 2-ethylhexyl ester |
| 16. Benzene, 1-ethyl-2-methyl- | 41. 1-Hexanol, 2-ethyl- |
| 17. Benzene, 1-ethyl-3-methyl- | 42. Tetramethylbutanedinitrile |
| 18. 1-Methoxy-2-propyl acetate | 43. Acetophenone |
| 19. Benzene, 1,3,5-trimethyl- | 44. Pentanedioic acid, dimethyl ester |
| 20. 2-Heptanone, 4,6-dimethyl- | 45. Benzenemethanol, dimethyl- |
| 21. Styrene | 46. Propanoic acid, 2-methyl-, 3-hydroxy-2,4,4-trimethylpentyl ester |
| 22. Benzene, 1-ethyl-2-methyl- | 47. 2,2,4-Trimethyl-1,3-pentanediol diisobutyrate |
| 23. Benzene, 1,3,5-trimethyl- | 48. 2H-Pyran-3(4H)-one, 6-ethenyldihydro-2,2,6-trimethyl- |
| 24. Propanoic acid, 2-hydroxy-, 2-methylpropyl ester | 49. Tetramethyl-1,4-benzenedimethanol |
| 25. 9-Eicosyne | 50. Dibutyl phthalate |

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
Shinjuku-ku, Tokyo
163-1130, Japan

Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
USA

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
5652AM, Eindhoven
The Netherlands

Phone: +31-40-254-9531
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
Far East International Plaza
No.317 Xianxia Road, Changning District
Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors