

Concentration Analysis of Volatile Constituents of Candy
- Comparative Effects of MonoTrap RGC18 TD Collection

MonoTrap RGC18 TD and HandyTD TD265 were used to screen for volatile components in candy. When collecting volatile compounds, the differences in the effects of collection methods were compared depending on the condition of the samples.

MonoTrap RGC18 TD uses a trapping agent to collect volatile compounds, the HandyTD TD265 is used to introduce those compounds into a GC by thermal desorption. In this application, the collection of volatile components from whole candy, was compared with a method of crushing samples prior to collection. The sensitivity of volatile compounds was greatly improved when the sample was crushed rather than sampled from the whole sample. The finer sample may have increased the surface area and increased the efficiency of volatilization.

Pretreatment procedure

Candy (4.5 g per piece of candy)

Place a pellet of apple candy in 20 mL vials

Collection (HS)
MonoTrap RGC18 TD × 1
(Cat.No. :1050-74201)

37°C 1 hour
collection



Sample 1
(untreated)



Sample 2
(crushed)

HandyTD/GC/MS




(Cat.No. : 2709-80000)

GC/MS Conditions


System : Thermal Desorption-GC/MS
(HandyTD TD265)
Column : InertCap Pure-WAX
0.25 mm I.D. × 60 m, df = 0.5 μm
Col. Cat. No. : 1010 - 68164
Col.Temp. : 40 °C(5 min) - 10 °C/min - 250 °C
Carrier Gas : He, 1 mL/min (constant flow)
GC Inlet : 250 °C Split 10:1
Detection : MS Scan (*m/z*30-350)

HandyTD Conditions

Desorb Temp. : Room temperature- 45 °C/s-
200 °C(5 min)
Pre Desorb Press. : 140 kPa



Analytical results of candy (crushed)



* Standard samples are not used for qualitative analysis.
Results from a library search.

1. Methyl alcohol
2. Ethyl alcohol
3. Propyl acetate
4. Valeraldehyde
5. Isobutyl methyl ketone
6. Ethyl butyrate
7. Isopropyl 2-methyl butanoate
8. Ethyl isovalerate
9. Butyl acetate
10. Hexanal
11. Isoamyl acetate
12. Valeraldehyde diethyl acetal
13. Methyl caproate

14. Isoamyl alcohol
15. 2-Hexenal
16. Ethyl caproate
17. Hexyl acetate
18. Isoamyl 2-methyl butyrate
19. 2-Hexenyl acetate
20. Ethyl lactate
21. Hexanol
22. Hexanal propylene glycol acetal
23. 2-Hexen-1-ol
24. Hexenal propylene glycol acetal
25. Linalool oxide
26. Hexenal propylene glycol acetal

27. Linalool oxide
28. Benzaldehyde
29. Diethyl malonate
30. Propylene Glycol
31. Hexyl caproate
32. Citronellyl propionate
33. Geranyl propionate
34. Caproic acid
35. Damascenone
36. Phenylethyl alcohol
37. 2-Hexenoic acid
38. Triacetin
39. Decalactone

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly.
We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
Shinjuku-ku, Tokyo
163-1130, Japan

Phone: +81-3-5323-6620
Fax: +81-3-5323-6621
Email: world@glsc.co.jp
Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
Torrance, CA 90503
USA

Phone: +1-310-265-4424
Fax: +1-310-265-4425
Email: info@glsciencesinc.com
Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
5652AM, Eindhoven
The Netherlands

Phone: +31-40-254-9531
Email: info@glsciences.eu
Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
Far East International Plaza
No.317 Xianxia Road, Changning District
Shanghai, China 200051

Phone: +86-21-62782272
Email: contact@glsciences.com.cn
Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors