



Analysis of Sodium Valproate

(Under the Condition of the draft for the Japanese Pharmacopoeia,
Sodium Valproate Extended-release Tablets B)

Data No. LB532-0812



Conditions

System : GL7700 HPLC system
 Column : InertSustain C18 (5 µm, 150 x 4.6 mm I.D.)
 Column Cat. No. : 5020-07345
 Eluent : A) CH₃CN
 B) 50 mM NaH₂PO₄ (pH 3.0, H₃PO₄)
 A/B = 50/50, v/v
 Flow rate : 1.0 mL/min
 Col. Temp. : 25 °C
 Detection : UV 210 nm (UV7750 UV Detector)
 Sample : Standard

Analyte:

1. Methyl *p*-hydroxybenzoate 4 mg/L
 2. Valproic acid 220 mg/L (2-a)
 or 800 mg/L (2-b)

Theoretical plates (2-a) : 9,617 (≥ 3,000)
 Symmetry factor (2-a) : 1.10 (≤ 2.0)
 RSD of the peak area of
 2-a (%) (n=6) : 0.23 (≤ 1.0)

Resolution (1, 2-b) : 16.4 (≥ 7)
 RSD of the relative peak area of
 2-b to 1 (%) (n=6) : 0.07 (≤ 1.0)