

Analysis of Phthalic Acid


Effect of Mobile Phase pH on Retention Time

When the pH of mobile phase is changed in the analysis of reverse phase mode on HPLC, the dissociation / non-dissociation state changes depending on the compound, so the retention time changes. Normally, the retention time is longer in non-dissociated state and shorter in dissociated state. This time, using Inertsil WP300 C18, we analyzed the basic compound saccharin and the acidic compound phthalic anhydride standard, and found that the elution order changed under acidic and neutral conditions. We also report that it was estimated that phthalic anhydride in the sample was partially changed to phthalic acid. (R.Takahashi)

Measurement Examples of Saccharin and Phthalic Anhydride

Acidic conditions

- : A) CH₃OH
 B) 50 mM KH₂PO₄ (0.1% H₃PO₄ in H₂O)
 (pH 2.6)




The retention of acidic compounds phthalic acid and phthalic anhydride varies greatly depending on the acidic and neutral conditions.

The retention of the basic compound saccharin remains almost unchanged.



Neutral condition

- : A) CH₃OH
 B) 50 mM K₂HPO₄ (0.1% H₃PO₄ in H₂O)
 (pH 6.7)



For the analysis conditions, we referred to the 18th revised Japanese Pharmacopoeia draft (acidic conditions) and the 18th revised Japanese Pharmacopoeia (neutral conditions).

Phosphoric Acid Buffering Action and Phthalic Acid Dissociation




Phosphoric acid has a higher buffering capacity near pH = 2.15, 7.20, 12.35. Even with the same phosphate buffer solution, the retention of acidic compounds is strong under acidic conditions, while the retention of basic compounds is strong under basic conditions.

HPLC condition


Column : Inertsil WP300 C18
 5 μ m, 150 \times 4.6 mm I.D.
 Flow rate : 1.0 mL/min
 Column Temp. : 20 $^{\circ}$ C
 Detection : UV 230 nm
 Inj. Vol. : 10 μ L
 Gradient Cond. : The solvent is described on the previous page.

Time(min)	A(vol%)	B(vol%)
0	10	90
7.0	10	90
8.0	95	5
11.0	95	5


Chemical Structures



Phthalic acid



Phthalic anhydride




Saccharin

Structures are created using Chemistry 4-D Draw which is provided by ChemInnovation Software, Inc.

Phthalic anhydride standard product

When the standard phthalic anhydride was measured, two peaks were seen, and it was confirmed that they were phthalic acid and phthalic anhydride, respectively.



Comparing phthalic acid and phthalic anhydride, there was a difference in the absorption curves around 270 nm to 280 nm.

GL Sciences disclaims any and all responsibility for any injury or damage which may be caused by this data directly or indirectly. We reserve the right to amend this information or data at any time and without any prior announcement.

GL Sciences Inc. Japan

22-1 Nishishinjuku 6-chome
 Shinjuku-ku, Tokyo
 163-1130, Japan

Phone: +81-3-5323-6620
 Fax: +81-3-5323-6621
 Email: world@glsc.co.jp
 Web: www.glsciences.com

GL Sciences Inc. USA

4733 Torrance Blvd. Suite 255
 Torrance, CA 90503
 USA

Phone: +1-310-265-4424
 Fax: +1-310-265-4425
 Email: info@glsciencesinc.com
 Web: www.glsciencesinc.com

GL Sciences B.V.

Dillenburgstraat 7C
 5652AM, Eindhoven
 The Netherlands

Phone: +31-40-254-9531
 Email: info@glsciences.eu
 Web: www.glsciences.eu

GL Sciences (Shanghai) Limited

Tower B, Room 2003
 Far East International Plaza
 No.317 Xianxia Road, Changning District
 Shanghai, China 200051

Phone: +86-21-62782272
 Email: contact@glsciences.com.cn
 Web: www.glsciences.com.cn



International Distributors

Visit our Website at www.glsciences.com/distributors