

InertSearch™ for LC

Inertsil® Applications

Establishing an ECD(Electro Chemical Detector)-HPLC System Using a Unique Diamond Electrode High-precision Quantitative Analysis of SAA (Sulfur Amino Acids)

Data No. LL001-0000

The poster was presented at pittcon 2009.

Joint presentation with AJINOMOTO Co., Inc., Pharmaceutical Research Lab Japan.

Author name : Junichi Isegawa 1, Akira Nakayama 1, Naoko Arashida 1,
Izumi Miyazaki 2, Takao Tamura 2
1 AJINOMOTO CO., INC. Pharmaceutical Research Lab. 2 GL Sciences Inc.

Data source : poster

Year : 2009

Conditions

Column : Inertsil ODS-3

Detection : ECD (ED703 pulse EC Detector, Diamond)

Sample : Nutrition solutions
blood plasma

Analyte : L-Cysteine
L-Cystine